


# OAPL TRAVEL SOCKS

Reduce The Risk Of DVT



OAPL Travel Socks are clinically proven to reduce the risk of Deep Vein Thrombosis (DVT) forming in your legs during travel. This is particularly so after long periods of inactivity which can allow DVT to form due to sluggish circulation. DVT can migrate to the lungs with fatal results.

## Design

Wearing OAPL Travel Socks can reduce the risk of DVT forming by improving circulation. They can also reduce oedema in the lower leg. OAPL Travel Socks feature graduated compression (14-18mm HG) and are designed with a closed toe for patient comfort.

## Indications

Prevention of DVT, aching and swelling legs, foot and leg oedemas.

## Contra-Indications

Severe arteriosclerosis, skin inflammation or ulcers, recent surgical wounds, gross oedema, cardiac failure or gross leg deformities. Please consult your physician.

OAPL Travel Socks are latex free!



SIZE	BLACK BEIGE		MEASUREMENTS	
	CODE		ANKLE CIRC	LEG HEIGHT
Small	41032	41036	Up to 21cm	Up to 40cm
Medium	41033	41037	22cm - 26cm	Up to 42cm
Large	41034	41038	26cm - 29cm	Up to 45cm
X-Large	41035	41039	Over 30cm	Up to 48cm